



Biohit Tip Compatibility with Other Pipette Brands

Optifit Tips, non-filtered

SafetySpace™ Filter Tips

Innovating for Health

		10µl	10µl ext	300µl	350µl	1000µl	1000µl WB	1200µl	5000µl	10ml	10µl	10µl ext*	20µl	20µl G*	120µl	200µl	300µl	500µl*	1000µl	1200µl		
Eppendorf																						
Reference	0.1-2.5µl	.	.									-	.									
	0.5-10µl									
	2-20µl (gray)			-									-									
	2-20µl (yellow)			.									.									
	10-100µl			.	.											.						
	50-200µl			.	.											.						
100-1000µl						
500-2500µl									.												.	
Research	0.1-2.5µl	.	.									-	.									
	0.5-10µl									
	2-20µl			.								.	.									
	10-100µl								
	20-200µl								
	100-1000µl						
	500-5000µl								.												.	
1000-10000µl								.												.		
Research 8-ch	0.5-10µl									
	30-300µl			-	.													.				
Research Plus	20-200µl						
8-ch	30-300µl				.												.	.				
Gilson																						
Ultra	0.2-2µl	*	*									*	*									
	1-10µl	*	*									*	*									
	2-20µl			.										-	.							
	20-100µl			.	.										.							
	20-200µl						
	200-1000µl						.	.													.	
Pipetman P/ Pipetman Neo	0.2-2µl									
	1-10µl									
	2-20µl			.									-	.								
	10/20-100µl							
	20/50-200µl							
100/200-1000µl							
Pipetman P	1000-5000µl								-												.	
	1-10ml									.											.	
Pipetman Neo 8-ch	20-200µl			
Brand																						
TransferpetteS	0.5-10µl									
	2-20µl			.									.									
	10-100µl			-						
	20-200µl						
	100-1000µl						
	0.5-5ml								-												.	
Transferpette 8-ch	30-300µl				
Thermo																						
Focus (short)	0.3-3µl									
Focus (short)	0.5-5µl									
Focus (short)	100-1000µl						.	.													.	
Focus	10ml									-											.	
Focus 8-ch	1-10µl									
Focus 12-ch	1-10µl									
Focus 8-ch	30-300µl						
Focus 12-ch	30-300µl						
Focus 12-ch	5-50µl			.											-							
Finnpipette Digital	0.2-2µl									
Finnpipette Digital	0.5-10µl									
Finnpipette Digital	10-100µl			.	.										.							
Finnpipette Digital	40-200µl						
Finnpipette Digital	100-1000µl						
Finnpipette Digital	200-1000µl						
Finnpipette Digital	1-5ml									-											.	
Finnpipette Digital	2-10ml									-											.	
Finnpipette Digital 8-ch	0.5-10µl									
Finnpipette Digital 8-ch	5-50µl			.											-							
Finnpipette Digital 8-ch	50-300µl						
Thermo F1	0.2-2µl									
Thermo F1	0.5-5µl									
Thermo F1	1-10µl									
Thermo F1	2-20µl			.									.									
Thermo F1	10-100µl						
Thermo F1	20-200µl						
Thermo F1	30-300µl						
Thermo F1	100-1000µl						
Thermo F1	0.5-5ml										-										.	
Thermo F1	1-10ml										-										.	
Thermo F1 8-ch	1-10µl									
Thermo F1 8-ch	5-50µl						
Thermo F1 8-ch	10-100µl						
Thermo F1 8-ch	30-300µl						
Thermo F2	0.2-2µl									
Thermo F2	0.5-5µl									
Thermo F2	1-10µl									
Thermo F2	2-20µl			.									.									
Thermo F2	10-100µl						
Thermo F2	20-200µl						
Thermo F2	100-1000µl						
Thermo F2	0.5-5ml										-										.	
Thermo F2	1-10ml										-										.	
Thermo F2 8-ch	1-10µl									
Thermo F2 8-ch	5-50µl						
Thermo F2 8-ch	10-100µl						
Thermo F2 8-ch	30-300µl						

The compatibility chart above is given for guidance. Biohit is not responsible for any technical modifications or changes made to the pipettes listed. Due to Biohit's continuing R&D effort, specifications may change without prior notice.

The specifications (=inaccuracy and imprecision values) are not tested.

The symbol indicates the fitting and the tip ejection of the Biohit tips as follows:
 • Fitting - No fitting * Tip ejection does not work

*: Filter Tips with regular air space
 G*: This 20µl tip is compatible with Gilson pipettes as shown on the chart, ordering code 783209.